MRAD 256 – Radiographic Technique 2

Lab Final Exam Marking Details

The timer (30 minutes) starts when the student either selects their patient or when they begin writing notes on the paper provided.

- Students may use tracking when lining up the tube and detector. Due to the large angles, it appears off. The student MUST tell the instructor they are using tracking before the image is checked.
- Students can also align the tube to the line on the IR.
- Reasonable collimation not to exceed ~2.5 cm 1 inch past required anatomy.
- The student being tested is responsible for making sure the patient is moving as requested. However if the patient is intentionally helping or hindering the student being tested, that student will receive a 0 on their own exam.
- At the instructor's discretion, an "emergency pause" on the timer can be used if the student is upset and needs.

to compose themselves. The student may not refer to the requisition or any notes they wrote during this time.

- Under the positioning category, the following errors also belong in the category with "patient incorrectly positioned."
    - o Incorrect tube angle
    - o Incorrect centering causing clipped anatomy